

Table 3 Comprison between the multiplex RT-PCR (M) and the conventional RT-PCR (C) for detection of PRRSV blood samples from swine artificial infected

Item	Sample No.	Days post infection									
		2 <sup>th</sup>	4 <sup>th</sup>	6 <sup>th</sup>	8 <sup>th</sup>	10 <sup>th</sup>	12 <sup>th</sup>	14 <sup>th</sup>	16 <sup>th</sup>	18 <sup>th</sup>	20 <sup>th</sup>
M	No.1	-	+	+	+	+	+	+	-	-	-
	No.2	-	-	+	+	+	+	+	+	-	-
	No.3	-	-	+	+	+	+	+	+	-	-
	Control	-	-	-	-	-	-	-	-	-	-
C	No.1	-	-	-	+	+	+	+	-	-	-
	No.2	-	-	-	+	+	+	+	-	-	-
	No.3	-	-	+	+	+	+	+	-	-	-
	Control	-	-	-	-	-	-	-	-	-	-

Note: +: Positive; -: Negative